

BASIS HI

Setting a New Industry Standard

ACRYLIC



- High Impact Formula
- Optimal Dimensional Stability
- Exceptional Flexural Strength
- 1 ½ hours Cure Cycle
- Japanese Technology


BASIS HI



www.ydp-usa.com
www.yamahachidentalusa.com
Phone: (877) 320-2909
info@yamahachidentalusa.com

yamahachi 
Dental Products USA

BASIS HI is a Brand distributed by Yamahachi Dental Products USA Inc.



BASIS HI MAIN FEATURES

- High Impact Formula
- Optimal Dimensional Stability
- Exceptional Flexural Strength
- 1 ½ hours Cure Cycle
- Japanese Technology



ACRYLIC BASIS HI TECHNICAL INSTRUCTIONS

CURING INSTRUCTIONS

1.	Powder / Liquid Ratio gm/ml 100:43 Immerse the flask in a container of tap water.
2.	Apply heat until boil. Let the resin completely cure for 30 – 40 minutes (Curing time starts when the water with the flask has started to boil).
3.	Cool the flask for about 30 minutes at room temperature.
4.	Recover denture after cooling completely.

PHYSICAL PROPERTIES

Flexural Strength, MPa	112.7
Flexural Modulus, MPa	2400
Vickers Hardness, Hv	19.5
Residual Monomer wt %	0.7
Shades:	O Pink, V Pink, LF Pink 

Shade Guides Available Upon Request.

<https://yamahachidentalusa.com/products/basis-hi-acrylic-resin>

Package information:

Powder: 1 Kilo and 3 Kilos
Liquid: 500 ml


BASIS HI



www.ydp-usa.com
www.yamahachidentalusa.com
Phone: (877) 320-2909
info@yamahachidentalusa.com

yamahachi 
Dental Products USA

BASIS HI is a Brand distributed by Yamahachi Dental Products USA Inc.