



SHINE 3D ML

The centerpiece of our Elements Zirconium range is the SHINE 3D PRO Multilayer disc. Offering exceptional translucency, from 42% to 49%, it is distinguished by a quality far superior to other discs.

SHINE 3D is designed with a unique shade layering and progressive hardness, ranging from 700 MPa in the incisal area to 1100 MPa at the cervical level.

Specially developed for monolithic anterior restorations, high-quality aesthetic veneers and natural-looking teeth, SHINE 3D is manufactured from rigorously selected raw materials. Available in shades A1 to BL4.

Description

Composition

ZrO₂ + HfO₂ + Y₂O₃: ≥99,0%

Technical Data

Flexural strength	700 MPa at the incisal to 1100 MPa at the cervical part
CTE	10,5
Density	6,07g/cm ³

A1

Item Nr.: 311-142112	12 mm
Item Nr.: 311-142114	14 mm
Item Nr.: 311-142116	16 mm
Item Nr.: 311-142118	18 mm
Item Nr.: 311-142120	20 mm
Item Nr.: 311-142122	22 mm
Item Nr.: 311-142125	25 mm
A2	
Item Nr.: 311-142212	12 mm
Item Nr.: 311-142214	14 mm
Item Nr.: 311-142216	16 mm
Item Nr.: 311-142218	18 mm
Item Nr.: 311-142220	20 mm
Item Nr.: 311-142222	22 mm
Item Nr.: 311-142225	25 mm
A3	
Item Nr.: 311-142312	12 mm

Item Nr.: 311-142314	14 mm
Item Nr.: 311-142316	16 mm
Item Nr.: 311-142318	18 mm
Item Nr.: 311-142320	20 mm
Item Nr.: 311-142322	22 mm
Item Nr.: 311-142325	25 mm
A3,5	
Item Nr.: 311-142412	12 mm
Item Nr.: 311-142414	14 mm
Item Nr.: 311-142416	16 mm
Item Nr.: 311-142418	18 mm
Item Nr.: 311-142420	20 mm
Item Nr.: 311-142422	22 mm
Item Nr.: 311-142425	25 mm
A4	
Item Nr.: 311-142512	12 mm
Item Nr.: 311-142514	14 mm

Item Nr.: 311-142516	16 mm
Item Nr.: 311-142518	18 mm
Item Nr.: 311-142520	20 mm
Item Nr.: 311-142522	22 mm
Item Nr.: 311-142525	25 mm
B1	
Item Nr.: 311-142612	12 mm
Item Nr.: 311-142614	14 mm
Item Nr.: 311-142616	16 mm
Item Nr.: 311-142618	18 mm
Item Nr.: 311-142620	20 mm
Item Nr.: 311-142622	22 mm
Item Nr.: 311-142625	25 mm
B2	
Item Nr.: 311-142712	12 mm
Item Nr.: 311-142714	14 mm
Item Nr.: 311-142716	16 mm

Item Nr.: 311-142718	18 mm
Item Nr.: 311-142720	20 mm
Item Nr.: 311-142722	22 mm
Item Nr.: 311-142725	25 mm
B3	
Item Nr.: 311-142812	12 mm
Item Nr.: 311-142814	14 mm
Item Nr.: 311-142816	16 mm
Item Nr.: 311-142818	18 mm
Item Nr.: 311-142820	20 mm
Item Nr.: 311-142822	22 mm
Item Nr.: 311-142825	25 mm
B4	
Item Nr.: 311-142912	12 mm
Item Nr.: 311-142914	14 mm
Item Nr.: 311-142916	16 mm
Item Nr.: 311-142918	18 mm

Item Nr.: 311-142920	20 mm
Item Nr.: 311-142922	22 mm
Item Nr.: 311-142925	25 mm
C1	
Item Nr.: 311-143012	12 mm
Item Nr.: 311-143014	14 mm
Item Nr.: 311-143016	16 mm
Item Nr.: 311-143018	18 mm
Item Nr.: 311-143020	20 mm
Item Nr.: 311-143022	22 mm
Item Nr.: 311-143025	25 mm
C2	
Item Nr.: 311-143112	12 mm
Item Nr.: 311-143114	14 mm
Item Nr.: 311-143116	16 mm
Item Nr.: 311-143118	18 mm
Item Nr.: 311-143120	20 mm

Item Nr.: 311-143122	22 mm
Item Nr.: 311-143125	25 mm
C3	
Item Nr.: 311-143212	12 mm
Item Nr.: 311-143214	14 mm
Item Nr.: 311-143216	16 mm
Item Nr.: 311-143218	18 mm
Item Nr.: 311-143220	20 mm
Item Nr.: 311-143222	22 mm
Item Nr.: 311-143225	25 mm
C4	
Item Nr.: 311-143312	12 mm
Item Nr.: 311-143314	14 mm
Item Nr.: 311-143316	16 mm
Item Nr.: 311-143318	18 mm
Item Nr.: 311-143320	20 mm
Item Nr.: 311-143322	22 mm

Item Nr.: 311-143325	25 mm
D2	
Item Nr.: 311-143412	12 mm
Item Nr.: 311-143414	14 mm
Item Nr.: 311-143416	16 mm
Item Nr.: 311-143418	18 mm
Item Nr.: 311-143420	20 mm
Item Nr.: 311-143422	22 mm
Item Nr.: 311-143425	25 mm
D3	
Item Nr.: 311-143512	12 mm
Item Nr.: 311-143514	14 mm
Item Nr.: 311-143516	16 mm
Item Nr.: 311-143518	18 mm
Item Nr.: 311-143520	20 mm
Item Nr.: 311-143522	22 mm
Item Nr.: 311-143525	25 mm

D4	
Item Nr.: 311-143612	12 mm
Item Nr.: 311-143614	14 mm
Item Nr.: 311-143616	16 mm
Item Nr.: 311-143618	18 mm
Item Nr.: 311-143620	20 mm
Item Nr.: 311-143622	22 mm
Item Nr.: 311-143625	25 mm
BL1	
Item Nr.: 311-143712	12 mm
Item Nr.: 311-143714	14 mm
Item Nr.: 311-143716	16 mm
Item Nr.: 311-143718	18 mm
Item Nr.: 311-143720	20 mm
Item Nr.: 311-143722	22 mm
Item Nr.: 311-143725	25 mm
BL2	
Item Nr.: 311-143812	12

	mm
Item Nr.: 311-143814	14
	mm
Item Nr.: 311-143816	16
	mm
Item Nr.: 311-143818	18
	mm
Item Nr.: 311-143820	20
	mm
Item Nr.: 311-143822	22
	mm
Item Nr.: 311-143825	25
	mm
BL3	
Item Nr.: 311-143912	12
	mm
Item Nr.: 311-143914	14
	mm
Item Nr.: 311-143916	16
	mm
Item Nr.: 311-143918	18
	mm
Item Nr.: 311-143920	20
	mm
Item Nr.: 311-143922	22
	mm
Item Nr.: 311-143925	25
	mm
BL4	
Item Nr.: 311-144012	12
	mm
Item Nr.: 311-144014	14

	mm
Item Nr.: 311-144016	16
	mm
Item Nr.: 311-144018	18
	mm
Item Nr.: 311-144020	20
	mm
Item Nr.: 311-144022	22
	mm
Item Nr.: 311-144025	25
	mm

