

LED CURE Master

[Technology](#)[Specifications](#)[Leaflet](#)

Specifications

Usage environment	Temperature 5-40°C, Humidity 10-95%RH, altitude max2000m
Operating Voltage and Frequency	AC230V±10% 50Hz
Power consumption	160VA
Number of course	4 Courses
Curing time	5-995 Sec.
Drying time	5-995 Sec. (Course button 4.)
Revolutions per minute	8.3r/min(50Hz)
External dimensions	210(W)x220(H)x223(D)mm
Dimensions of chamber	135(W)x77(H)x135(D)mm
Weight	5.4kg
Emission wavelength	375nm ~ 495nm
Quantity of LED lamps	Power LED lamp 40pcs
Accessories	Tray: 1 pc, Big Pin: 5 pcs, Small Pin: 5 pcs

